


SKATE Plate

Cervical Plate Kit






Fusion with more Stability

Stabilization plate preventing the migration of the cage


SKATE Cervical Plate Kit

The SKATE Plate is a safe and effective device.

The plate is designed to avoid the complications related to use of cervical intrvertebral cage: risk of cage migration and inability to maintain the spinal stability.



Anatomic Fit



■ Anatomical Fixation

The plate can be bended, twisted to match anatomical curves without any strength against bone.

■ Simplicity

- 2-in-1 system
- Easy and quick insertion
- Secured fixation of the spinal segment

■ Multi-level Implantation

The SKATE Plate allows to combine several plate-cages in multi-segment arthrodesis.

■ Locking

Special plate hole design which prevents screw from pulling out.

Specification

REF	Width (mm)	Height (mm)
2100-0840P	14	4
2100-0940P	14	5
2100-1040P	14	6
2100-1140P	14	7
2100-1240P	14	8



INSTRUMENTS

Cat. Number Description

218-0624 Anatomic Trial No.4



218-0625 Anatomic Trial No.5



218-0626 Anatomic Trial No.6



218-0627 Anatomic Trial No.7




218-0628 Anatomic Trial No.8



218-0201 3.5 Awl




218-0525 2.5 Hex screw driver



218-0807 Cervical Cage Holder



218-0810 SKate Plate Holder



218-0003 Instrument Tray

