


# SenCho

## Sacroiliac System

SenCho Sacroiliac System consists of two unique diameter cannulated implants offered in various lengths, in full and partially threaded designs to accommodate different patient anatomy.

The system is intended for sacroiliac joint fusion for conditions including sacroiliac joint disruptions and degenerative sacroiliitis.




### Features

#### SI Screw


##### A Dual-thread Pattern

Enhance compression effect and stability




##### B Rough Surface

Increase the frictional Force



##### C Fenestration Design

Promote bone fusion



#### Compression Screw


##### A Compression Screw and Washer Combination

Increase the compression force between ilium and sacrum to promote SI joint combination



##### B Free Washer Design

- Fit the ilium
- Improve load distribution




# Specification

## SI Screw (Titanium)



REF	Dia. (D) (mm)	Length (L) (mm)	REF	Dia. (D) (mm)	Length (L) (mm)
21060010003000	10	30	21060011503000	11.5	30
21060010003500	10	35	21060011503500	11.5	35
21060010004000	10	40	21060011504000	11.5	40
21060010004500	10	45	20160011504500	11.5	45
21060010005000	10	50	21060011505000	11.5	50
21060010005500	10	55	21060011505500	11.5	55
21060010006000	10	60	21060011506000	11.5	60
21060010006500	10	65	21060011506500	11.5	65
21060010007000	10	70	21060011507000	11.5	70
21060010503000	10.5	30	21060012003000	12	30
21060010503500	10.5	35	21060012003500	12	35
21060010504000	10.5	40	21060012004000	12	40
20160010504500	10.5	45	20160012004500	12	45
21060010505000	10.5	50	21060012005000	12	50
21060010505500	10.5	55	21060012005500	12	55
21060010506000	10.5	60	21060012006000	12	60
21060010506500	10.5	65	21060012006500	12	65
21060010507000	10.5	70	21060012007000	12	70
21060011003000	11	30	21060012503000	12.5	30
21060011003500	11	35	21060012503500	12.5	35
21060011004000	11	40	21060012504000	12.5	40
20160011004500	11	45	20160012504500	12.5	45
21060011005000	11	50	21060012505000	12.5	50
21060011005500	11	55	21060012505500	12.5	55
21060011006000	11	60	21060012506000	12.5	60
21060011006500	11	65	21060012506500	12.5	65
21060011007000	11	70	21060012507000	12.5	70

## Compression Screw (Titanium, Washer Included)



REF	Dia. (D) (mm)	Length (L) (mm)
21060107003001	7	30
21060107003501	7	35
21060107004001	7	40
21060107004501	7	45
21060107005001	7	50
21060107005501	7	55
21060107006001	7	60
21060107006501	7	65
21060107007001	7	70